Lab Use Only				
Res Lab #:				
Date Rec'd:				
Initials:				

DNA

National Inherited Bleeding Disorder Genotyping Laboratory



Department of Pathology and Molecular Medicine Queen's University, Kingston, Ontario

Hemophilia A and B Genotype Testing Requisition

Patient Name: Surname, First Name) Sex: Male D Female D					
DOB://	Unique Identifier:	Health card #, Hosp	Dital # CBDR #:		
Date of specimen collection:/// Phlebotomist:					
Referring Clinic:	Report to	:	Fax #: _		
Test Requested:	Hemophilia A 🛛		Hemophilia B 🛛		
Coagulation Factor Level:	Factor VIII	U/mL	Factor IX	U/mL	
Inhibitor: Yes 🗆 No 🗆	Inhibitor Titre	:	B.U.		
Has intron 22 inversion testing been done? Yes D No D					
Information Requested:	 Confirmation of dia Carrier status Prenatal diagnosis 	0			
Pregnant? Yes D No D					
Have samples from this family been sent to this lab before? Yes \Box No \Box					
If Yes, specify					
Relationship to this patient					
Sample Requirements: Ship to:					
6 cc whole blood		Attn: Gina Jones/Samira Kheitan			
EDTA (lavender top) <u>or</u>		Department of Pathology and Molecular Medicine Queen's University, Richardson Laboratory, Room 201 88 Stuart St., Kingston, Ontario K7L 3N6			
ACD (yellow top) <u>or</u>					

ieen's University, Richardson Laboratory, Room 88 Stuart St., Kingston, Ontario K7L 3N6 Tel: 613-533-3187 FAX: 343-344-2733 Email: NIBDGL@ queensu.ca

National Inherited Bleeding Disorder Genotyping Lab Sample Collection Instructions

Requisitions:

- 1. Samples must be accompanied by a completed requisition form.
- 2. Submitted patients must have a documented rationale for testing:
 - For affected hemophilia A or B patients, a clotting factor activity level must be provided.
 - For VWD patients, VWF:Ag, VWF activity (specify test used), FVIII:C, and multimer information must be provided.
 - For carrier testing, a documented family history or a coagulation factor level must be provided.
 - For carrier testing, information on the family variant if available must be provided.
 - For prenatal testing, information on the family variant if available must be provided.
- 3. Incomplete requisition forms will result in delayed sample testing.

Sample Collection and Shipment:

- 1. Samples acceptable for testing:
 - Venous whole blood (minimum 6 cc) collected into EDTA (lavender top) or ACD (yellow top) evacuated tubes
 - Expired tubes should not be used.
 - If blood is being drawn from an intravenous line for laboratory testing, two times the deadspace volume should be discarded.
 - When drawing blood specimens for several examinations during a single venipuncture, the "order of draw" shall be: (1) blood culture; (2) coagulation specimens; (3) serum tube with or without clot activator or gel; (4) heparin; (5) EDTA; (6) glycolytic inhibitor.
 - DNA (minimum of 15 µg at 150 ng/µL); for some cases smaller samples are acceptable.
 - Patient ID must be verified and samples labelled using two unique identifiers.
 - Samples shall be collected using routine practices/standard precautions.
 - Materials for sample collection shall be safely disposed of according to institutional protocols.
- 2. Ship packages on a Monday, Tuesday, Wednesday, Thursday as follows:
 - Shipping temperature: DNA (ambient), whole blood (cold packs), whole blood frozen (dry ice).
 - Place the samples in sealable plastic bags with absorbent material.
 - Include completed requisition and consent form.
 - Ship overnight. Contact courier for complete shipping instructions.
- 3. Attach the following labels to the outside of the box/package:
 - Dry ice label (if applicable).
 - Return address label (including the contact name and telephone number).
- 4. Ship to: Attn: Gina Jones/Samira Kheitan
 - Department of Pathology and Molecular Medicine Queen's University, Richardson Laboratory, Room 201 88 Stuart Street, Kingston, Ontario K7L 3N6
 - Tel: 613-533-3187 FAX: 343-344-2733 NIBDGL@queensu.ca

Results:

- Turnaround time is approximately 3 months from the time of sample submission but may be longer for rare genes or for carrier testing if familial variant is unknown.
- If the family-specific variant is known, urgent reporting can be completed within several weeks. Please indicate this on the requisition form.
- For prenatal testing, maternal cell contamination studies are the responsibility of the referring clinic.
- Results will be faxed to the referring clinician or designate listed on the requisition.
- Results will also be entered into the Canadian Bleeding Disorders Registry (CBDR) when possible.